Right out of the box, it lives up to its name.


INTRODUCING THE 3M™ LITTMANN® CLASSIC III™ STETHOSCOPE. THE LATEST VERSION OF THE STETHOSCOPE THAT HELPS MILLIONS OF MEDICAL PROFESSIONALS ACHIEVE THEIR BEST.

The Classic III stethoscope brings new design, materials, and technology to the Littmann stethoscope series that's been used and trusted by medical professionals for decades. With a two-sided chestpiece, dual tunable diaphragms, improved tubing, and much more, the Classic III stethoscope combines the best of new and traditional.

Made in the USA, the Classic III stethoscope is a durable, dependable, clinical tool that offers high acoustic sensitivity for performing general physical assessments. Using a Classic III stethoscope demonstrates what's important to you as a medical professional: quality, achievement, and providing the best possible care. Every patient. Every day.



It's an investment you make in yourself.





3M Health Care


3M Infection Prevention Division

1-800-228-3957 | www.littmann.com

2510 Conway Ave., St. Paul, MN 55144 USA

3M, Littmann, the Littmann Quality seal, the L, and

Please recycle. Printed in USA. © 3M 2014. All rights reserved. Made in USA.



3M™ Littmann® Classic III™

STETHOSCOPE

THE CLASSIC III STETHOSCOPE OFFERS HIGH ACOUSTIC SENSITIVITY FOR EXCEPTIONAL PERFORMANCE WHEN DOING GENERAL PHYSICAL ASSESSMENTS.



Tunable diaphragms on both the adult and pediatric sides of the chestpiece. The pediatric side is useful for small or thin patients, around bandages, and for carotid assessment.



Open bell stays clear of dirt and debris by covering it with the small diaphragm.



No natural rubber latex or phthalate plasticizers are used in the tubing or any other component.



Next-generation tubing

lasts longer with improved resistance to skin oils and alcohol; less likely to pick up stains.



Single-piece diaphragm is

easy to attach and easy to clean because the surface is smooth without crevices.



Open bell versatility

The paediatric side of the Classic III stethoscope converts to a traditional open bell. Simply remove the single-piece diaphragm and replace it with a non-chill rim.

