

3M™ Littmann® Stethoscopes Product Comparison Chart

An investment you make in yourself.

Littmann stethoscopes are made for professionals who are dedicated to achieving the best patient outcomes. They offer superb acoustic performance, consistently high quality, and an outstanding value.



Lightweight		Classic				Cardiology			Electronic	
<p>6</p> <p>Acoustic Performance</p> <p>Ideal for taking blood pressure readings and making limited patient assessments.</p>		<p>7-8</p> <p>Acoustic Performance</p> <p>Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.</p>				<p>9-10</p> <p>Acoustic Performance</p> <p>Designed to meet the needs of demanding clinicians. Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation.</p>			<p>10+</p> <p>Acoustic Performance</p> <p>High-performance auscultation. Hear more of what you need to hear with amplification, ambient noise reduction, and optional connectivity.</p>	
		<p>NOW AVAILABLE</p> <p>The Next Generation</p>								
Lightweight II S.E. Stethoscope	Select Stethoscope	Classic II Pediatric and Infant Stethoscopes	Classic II S.E. Stethoscope	Classic III™ Stethoscope	Master Classic II™ Stethoscope	Cardiology III™ Stethoscope	Cardiology S.T.C. Stethoscope	Master Cardiology™ Stethoscope	Electronic Stethoscope Model 3100	Electronic Stethoscope Model 3200
Double sided chestpiece 118 grams Tunable diaphragm 2 year warranty	Single sided chestpiece 120 grams Tunable diaphragm 2 year warranty	Double sided chestpiece 105, 95 grams 3 year warranty	Double sided chestpiece 135 grams Tunable diaphragm 3 year warranty	Double sided chestpiece 150 grams Two tunable diaphragms Convertible open bell 5 year warranty	Single sided chestpiece 160 grams Tunable diaphragm 3 year warranty	Double sided chestpiece 175, 180 grams Available in 22" / 56 cm and 27" / 69 cm Two tunable diaphragms Dual-lumen tubing 5 year warranty	Single sided chestpiece with soft touch insert 180 grams Tunable diaphragm Dual-lumen tubing 6 year warranty	Single sided chestpiece 175, 185 grams Available in 22" / 56 cm and 27" / 69 cm Tunable diaphragm Dual-lumen tubing 7 year warranty	185 grams 2 year warranty	Bluetooth® wireless technology Recording and playback 185 grams 2 year warranty

Tunable Technology

3M invented tunable diaphragm technology to let you hear different frequency sounds without turning over and repositioning the chestpiece.



Light contact:

Hold the tunable diaphragm side of the chestpiece with light pressure to hear low frequency sounds.



Firm pressure:

Press the chestpiece a little more firmly to hear higher frequency sounds.

Quality you can trust

All Littmann mechanical stethoscopes are designed and manufactured in the United States according to strict process and quality controls. The result: consistent quality. That means any Littmann mechanical stethoscope you pick up will provide a consistently high level of performance.